

Catalog #201192B (500 µL)

Cell-ID™ Intercalator-Ir 500 µM



WARNING! CHEMICAL HAZARD. Before handling any chemicals, refer to the Material Safety Data Sheet (MSDS) provided by the manufacturer, and observe all relevant precautions.



NOTICE: HIGH CONCENTRATION. Cell-ID Intercalator-Ir 500 µM is a highly concentrated metal intercalator solution and must be diluted in accordance with this protocol to avoid early failure of the detector.

Description

Cell-ID Intercalator-Ir is a cationic nucleic acid intercalator that contains natural abundance Iridium (^{191}Ir and ^{193}Ir) and is used for identifying nucleated cells in CyTOF® analysis. When cells are stained with Intercalator-Ir, it will bind to cellular nucleic acid, and detection of both stable isotopes will enable identification of nucleated cells. It is a live cell membrane-impermeable dye and therefore requires cells to be fixed and/or permeabilized before staining.



Note: While dilutions of the 500 µM stock solution are suggested in the protocols below, the concentration can be titrated for individual cell types and experiments for optimal Cell-ID Intercalator staining. It is suggested not to exceed 1 µM intercalator concentration in the staining solution.

Staining Protocol A

- 1 Before intercalating, cells must be fixed.
 - If fixed with methanol, wash cells with PBS (without Ca^{2+} or Mg^{2+}) before proceeding.
 - Cells may be used directly if fixed with formaldehyde (3.7%, 30min, RT).
- 2 Dilute Cell-ID Intercalator-Ir 1:2000 with PBS (without Ca^{2+} or Mg^{2+}).
- 3 Use 0.5mL of working solution per 1×10^6 cells/ tube.
- 4 Incubate 15-20 mins at room temperature.
- 5 Wash cells with 2 mL PBS (without Ca^{2+} or Mg^{2+}) per tube. Repeat once.

Staining Protocol B (for use with the MaxPar[®] Cell Surface Staining Protocol)

- 1 After cell staining is complete, prepare 1 ml of cell intercalation solution for each sample by diluting Cell-ID Intercalator-Ir 1:4000 into MaxPar[®] Fix and Perm Buffer (Fluidigm Cat. 201067) and mix by vortexing.
- 2 Add 1 ml of the intercalation solution prepared in step 1 to each tube and gently vortex. Incubate for 1 hour at room temperature or leave overnight at 4 °C.



Note: Cells can be left at 4 °C in the intercalation solution up to 48 hours.

- 3 Wash cells by adding 2 ml of MaxPar[®] Cell Staining Buffer (Fluidigm Cat. 201068), centrifuge and discard supernatant by aspiration.
- 4 Repeat for a total of two washes with MaxPar Cell Staining Buffer.
- 5 Wash cells with 2 ml of MaxPar[®] Water (Fluidigm Cat. 201069), centrifuge and discard supernatant by aspiration.
- 6 Leave cells pelleted until ready to run on CyTOF. Immediately prior to CyTOF data acquisition, adjust cell concentration to $2.5\text{-}5 \times 10^5/\text{ml}$ with MaxPar Water and filter cells into cell strainer cap tubes.
- 7 Acquire data on CyTOF.

Technical Support

Phone	Email
In the United States: 1.866.FLUIDLINE (1.866.358.4354) Outside the United States: 1.650.266.6100	U.S. and countries not in Europe or Asia: techsupport@fluidigm.com Europe: techsupporteurope@fluidigm.com Asia: techsupportasia@fluidigm.com

For Research Use Only. Not for use in diagnostic procedures.

Fluidigm, the Fluidigm logo, Cell-ID, CyTOF, and MaxPar are trademarks or registered trademarks of Fluidigm Corporation in the U.S. and/or other countries. All other trademarks are the property of their respective owners. © 2014 Fluidigm Corporation. All rights reserved.

Corporate Headquarters

7000 Shoreline Court, Suite 100
South San Francisco, CA 94080 USA
Toll-free: 1.866.FLUIDLINE | Fax: 650.871.7152
www.fluidigm.com

Sales

North America | +1 650.266.6170 | info-us@fluidigm.com
Europe/EMEA | +33 1 60 92 42 40 | info-europe@fluidigm.com
Japan | +81 3 3555 2150 | info-japan@fluidigm.com
Asia | +1 650.266.6000 | info-asia@fluidigm.com

