

Hyperion Tissue Imager

SITE REQUIREMENTS GUIDE



For Research Use Only. Not for use in diagnostic procedures.

Information in this publication is subject to change without notice. It is Fluidigm policy to improve products as new techniques and components become available. Therefore, Fluidigm reserves the right to change specifications at any time. Every effort has been made to avoid errors in the text, diagrams, illustrations, figures, and screen captures. However, Fluidigm assumes no responsibility for any errors or omissions. In no event shall Fluidigm be liable for any damages in connection with or arising from the use of this publication.

Patent and Limited License Information

Fluidigm products are covered by issued and pending patents in the United States and other countries. Patent and limited license information is available at fluidigm.com/legalnotices

Trademarks

Fluidigm, the Fluidigm logo, CyTOF, Helios, and Hyperion are trademarks or registered trademarks of Fluidigm Corporation in the United States, and/or other countries. All other trademarks are the sole property of their respective owners.

For EU's WEEE directive information, go to fluidigm.com/compliance

Complies with 21 CFR 1040.10 and 1040.11 except for deviations pursuant to Laser Notice No. 50, dated June 24, 2007.

© 2017 Fluidigm Corporation. All rights reserved. 09/2017

PN 400300

IMC-SRG-02 rev 5

For technical support visit fluidigm.com/support

EMAIL

North America support.northamerica@fluidigm.com

Europe support.europe@fluidigm.com

Asia techsupportasia@fluidigm.com

Latin America techsupportlatam@fluidigm.com

All other countries techsupport@fluidigm.com

PHONE

United States (toll-free) +1 866 359 4354

Europe +33 1 60 92 42 40

Japan +81 3 3662 2150

China (excluding Hong Kong) +86 21 3255 8368

All other countries +1 650 266 6100

Contents

About This Guide	4
Safety Alert Conventions	4
Safety Alerts for Chemicals	4
Safety Alerts for Instruments	4
Introduction and Specifications	6
Introduction	6
Cover Safety Interlock	7
Solid-State Laser Specifications	7
Hyperion Tissue Imager Specifications	8
Requirements	9
Space Requirements	9
Laser Ablation System Packaging	10
Crate Information	10
Environmental Requirements	11
Electrical Requirements	11
Gas Requirements	13
Argon Gas	13
Helium Gas	14
Exhaust Requirements	16
Additional Site Recommendations	16
Institutional Network connectivity	16
Internet connectivity	16
Telephone access	16
Appendix A: Safety	17
General Safety	17
Instrument Safety	17
Instrument Symbols	18
Chemical Safety	19

About This Guide



CAUTION ABBREVIATED SAFETY ALERTS. Hazard symbols and hazard types specified in procedures may be abbreviated in this document. For complete safety information, see Appendix A.

Safety Alert Conventions

This guide uses specific conventions for presenting information that may require your attention. Refer to the following safety alert conventions.

Safety Alerts for Chemicals

Fluidigm follows the United Nations Globally Harmonized System of Classification and Labelling of Chemicals (GHS) for communicating chemical hazard information. GHS provides a common means of classifying chemical hazards and a standardized approach to chemical label elements and safety data sheets (SDSs). Key elements include:

- Pictograms that consist of a symbol on a white background within a red diamond-shaped frame. Refer to the individual SDS for the applicable pictograms and warnings pertaining to the chemicals being used.



- Signal words that alert the user to a potential hazard and indicate the severity level. The signal words used for chemical hazards under GHS:

DANGER Indicates more severe hazards.

WARNING Indicates less severe hazards.

Safety Alerts for Instruments

For hazards associated with instruments, this guide uses the following indicators:

- Pictograms that consist of a symbol on a white background within a black triangle-shaped frame.



- Signal words that alert the user to a potential hazard and indicate the severity level. The signal words used for instrument hazards:

DANGER Indicates an imminent hazard that will result in severe injury or death if not avoided.

WARNING Indicates a potentially hazardous situation that could result in serious injury or death.

CAUTION Indicates a potentially hazardous situation that could result in minor or moderate personal injury.

IMPORTANT Indicates information necessary for proper use of products or successful outcome of experiments.

Safety Data Sheets

Read and understand the SDSs before handling chemicals. To obtain SDSs for chemicals ordered from Fluidigm Corporation, either alone or as part of this system, go to fluidigm.com/sds and search for the SDS using either the product name or the part number.

Some chemicals referred to in this user guide may not have been provided with your system. Obtain the SDSs for chemicals provided by other manufacturers from those manufacturers.

Introduction and Specifications



WARNING Only Fluidigm trained personnel or Fluidigm Field Service Engineers may access the instrument when the doors are open or the covers have been removed.

Introduction

The Hyperion™ Imaging System is a high-parameter imaging system capable of analyzing up to 50 individual parameters at single-cell resolution in tissue sections. The Hyperion Imaging System comprises a Hyperion Tissue Imager, a Helios and the Hyperion software for data acquisition and visualization. Do not attempt to operate the instrument before reading this Site Requirements Guide for all the safety information.

The Hyperion Tissue Imager is classified as a Class 1 laser product containing a Class 4 laser encased in the optical enclosure. This is based on the maximum output power available for the intended use. Therefore the laser is exempt from user control measures. Your facility may have its own guidelines and recommendations on the safe use of lasers.




Figure 1: Instrument label.

Cover Safety Interlock

The two doors on the top cover of the Hyperion Tissue Imager are equipped with safety interlocks that are triggered by contact closure with the frame. If doors of the top cover are open, the interlock will be engaged to prevent the laser from firing.



WARNING The cover of the Hyperion Tissue Imager functions as a safety interlock. When the cover is off, the AC connection to the laser will be disengaged.

Solid-State Laser Specifications

The following are specifications of the solid state laser installed in the Hyperion Tissue Imager.

Table 1: Specifications for the solid state laser

Performance Specifications	
Laser Type	Nd:YAG
Wavelength	213 nm
Energy Output	$\leq 3 \mu\text{J}$
Optimal Repetition rate	100, 200 Hz

Hyperion Tissue Imager Specifications

Table 2: Operation specifications for the Hyperion Tissue Imager

Operation Specifications	
Power	200VA
Operating voltage	100-240 VAC
Operating frequency	50/60 Hz
Power dissipation, full load	< 75 W
Recommended Operating Temperature	18–25 °C
Optimal Temperature	21 °C ± 1 °C

Table 3: Instrument dimensions for the Hyperion Tissue Imager

Instrument Dimensions	
Height	134 (53) cm (in)
Operating Height	132 (51) cm (in)
Width	56 (22) cm (in)
Depth	56 (22) cm (in)
Weight	170 (375) kg (lbs)


Table 4: Input and output connections for the Hyperion Tissue Imager

Input and Output Connections	
Input	AC Power, Argon Gas, Helium Gas, Sample Slide
Output	Ablated Sample, USB 2.0 communication cable, USB 3.0 camera, Trigger (Data acquisition card)

Requirements

Space Requirements

The recommended general footprint for installation of the Hyperion Tissue Imager attached to the Helios system is illustrated below. Ensure there is adequate space for the Hyperion Tissue Imager instrument, the Helios instrument, computer, and accessories. The frame of the Hyperion Tissue Imager is 56 cm wide and 56 cm deep and should be placed directly in front of the Helios instrument. A Fluidigm field service engineer will be on hand to unpack and install the Hyperion Tissue Imager, the Helios instrument, and its accessories.



Dimensions Units in cm/inches

Figure 2: Footprint of the Hyperion Tissue Imager, Helios, and accessories.

Laser Ablation System Packaging

The Hyperion Tissue Imager is shipped with the following items:

- Sample coupling assemblies
- Calibration/tuning slides
- Monitor
- Power and communication cables
- Gas connection lines

Crate Information

The Hyperion Tissue Imager is shipped in one fully packaged crate and two boxes. A standard pump truck with minimum rating for 1,600 lb to 1 ton is recommended for moving the crates and box if necessary. Once you have received the shipment, store in a dry place not exposed to weather until the scheduled installation date. Table 4 provides the contents of the instrument crate and boxes.

Table 5: Components of the shipping crate of the Hyperion Tissue Imager.

Crate	Components						
Crate 1	Hyperion Tissue Imager						
Box 1	Monitor						
Box 2	Auxiliary power cable	Signal communication cables	Argon/Helium gas connection lines	Data acquisition card	Sample coupling assemblies	Tuning slide	Injector

Table 6: Crate dimensions of the Hyperion Tissue Imager.

Component	Depth cm (in)	Width cm (in)	Height cm (in)	Weight kg (lbs)
Crate Dimensions	135 (53)	74 (29)	155 (61)	347 (765)

Environmental Requirements

The laboratory where the Hyperion Tissue Imager is installed should meet the following environmental requirements.

Table 7. Environmental Requirements for the Hyperion Tissue Imager

Environmental Requirements	
Recommended Operating Temperature	18–25 °C
Optimal Temperature	21 °C ± 1 °C
Humidity	10–85 %
Vibration Environment	VC-A

IMPORTANT The Hyperion Imaging System should perform according to specifications in vibration environments at or above VC-A.

- Free of smoke, dust, and corrosive fumes
- Out of direct sunlight and away from heating system

Electrical Requirements

Be aware that high voltages are present in the Hyperion Tissue Imager control box once the AC power is switched on.

Table 8: Electrical specifications for the Hyperion Tissue Imager and accessories

Electrical Specification	
Power Supply	24 VDC (120 V max)
Operating Voltage	100–240 V AC
Peak Current (per circuit)	200 VA
Operating Frequency	50/60 Hz
Operating Mode	Single phase operation
Power Dissipation, Full Load	< 75 W

Power Cord Requirements

The power cord set supplied with the instrument will meet the requirements for your country. Power to the Hyperion Tissue Imager is to be delivered from one 2 A single-phase 100–240 V AC, 50/60 Hz auxiliary electrical branch circuit. The equipment is designed for connection to the 2 Amp Auxiliary outlet on the electrical panel on the right side of the Helios instrument.



WARNING Do not disable the grounding feature of this instrument. Always connect to the three prong cable. Do not replace detachable mains supply cords by inadequately rated cords.

Table 9: Plug information for the Hyperion Tissue Imager and accessories

Accessories	Voltage (V AC)	Current (A)
Hyperion Tissue Imager	100–240	<1



WARNING Be aware that high voltages are present in the laser box once AC power is switched on.

Gas Requirements

Argon Gas

The Hyperion Tissue Imager instrument has been designed to utilize argon and helium as the sample carrier gases. Argon is supplied by the make-up gas line on the Helios instrument.

Table 10: Argon Gas Specifications

Argon Gas Specification	
Purity	≥99.996%
Supply Pressure	85 ±5 psi
Flow Rate	20.5 L/min

IMPORTANT An average argon gas cylinder with a pressure of 2600 psi would last approximately 7 hours of continuous use of the Hyperion Imaging System. A standard 280 L liquid argon dewar would last approximately 140 hours of continuous system use.

NOTE We recommend installing a 350 psi pressure relief valve with the liquid argon Dewar.

Helium Gas

Research grade helium gas is recommended for use with the Hyperion Tissue Imager instrument.

Table 11: Helium Gas Specifications

Helium Gas Specification	
Purity	≥99.999% (Grade 5.0)
Supply Pressure	30 ±5 psi
Flow Rate	0.25 L/min



WARNING The laboratory should be able to accommodate a full-size compressed gas tank or have other provisions for supplying helium to the laser. The recommended helium supply pressure is approximately 30 ±5 psi. Ensure that all gas connections are tight prior to powering up the laser and the software.



WARNING We recommend installing an oxygen sensor in the room where the operator and gas storage are located.

NOTE The helium regulator must be capable of handling 1 SLPM constant flow, providing a pressure range of 0–50 psi and a ¼ inch Swagelok® tube adapter. A dual stage model regulator is recommended for more constant pressure.

Safe Handling of Gas Cylinders

The permanent installation of gas supplies is the responsibility of the user and should conform to local safety and building codes. Observe the following safety precautions when handling argon gas cylinders:

- Fasten all gas cylinders securely to an immovable bulkhead or a permanent wall.
- When gas cylinders are stored in confined areas, ventilation should be adequate to prevent dangerous accumulations. Move or store gas cylinders only in a vertical position with the valve cap in place.
- Locate gas cylinders away from heat or ignition sources, including heat lamps. Cylinders have a pressure relief device that will release the contents of the cylinder if the temperature exceeds 52 °C (125 °F).

- When storing cylinders external to a building, the cylinders should be stored so that they are protected against temperature extremes (including direct rays of the sun) and stored above ground on a suitable floor.
- Gas cylinders should be clearly marked to identify the contents and status (e.g., full, empty).
- Do not attempt to refill gas cylinders.
- Use only approved mechanical regulators and hose connectors. Left-hand thread fittings are used for fuel gas tank connections whereas right-hand fittings are used for oxidant and support gas connections.
- Arrange gas hoses away from foot traffic to avoid damage.
- Perform periodic gas leak tests by applying a soap solution to all joints and seals.

Regulators

Contact your gas supplier for recommendations on the Helium regulator that is compatible with the cylinders they provide.

IMPORTANT The regulator must be capable of handling 1 SLPM constant flow, providing 30 psi, low pressure output to a ¼ inch swagelok (or compatible) female fitting. A dual stage model regulator is recommended for more constant pressure.

Exhaust Requirements

The Hyperion Tissue Imager is designed to allow minimal venting of any gases from the left side of the instrument. There are no specific requirements for exhaust.

The information provided above covers site preparation requirements for the Hyperion Tissue Imager only. For complete information on the Helios component of the Hyperion Imaging System, please refer to:

- Helios, A CyTOF System User Guide PN 400100

Additional Site Recommendations

Institutional Network connectivity

Enabling network connectivity to the instrument workstation allows for more efficient transfer of instrument data through the Ethernet network port of the instrument workstation. Contact your site IT Department.

IMPORTANT Data should not be written directly to a Network Drive during sample acquisition and must be written directly to the instrument workstation. Data may only be transferred/shared via network connection upon completion of sample acquisition.

Internet connectivity

It is highly recommended that you have the Hyperion Imaging System computer on a network with internet connection so that the Fluidigm Support team can better assist you and your team and to facilitate remote diagnostics of the instrument. Please work with your IT department to ensure that internet connectivity can occur prior to installation. Access to the internet can greatly reduce instrument down time and unanticipated field service/repair costs. Under no circumstance will Fluidigm be able to access the instrument, or associated data stored on the workstation, without the end user's express permission.

Telephone access

If internet connection is unavailable and mobile telephone reception is poor, it is recommended that you provide a landline phone in the immediate vicinity of the instrument workstation.

Appendix A: Safety

General Safety

In addition to your site-specific safety requirements, Fluidigm recommends the following general safety guidelines in all laboratory and manufacturing areas:

- Use personal protective equipment (PPE): safety glasses, fully enclosed shoes, lab coats, and gloves.
- Know the locations of all safety equipment (fire extinguishers, spill kits, eyewashes/showers, first-aid kits, safety data sheets, etc.), emergency exit locations, and emergency/injury reporting procedures.
- Do not eat, drink, or smoke in lab areas.
- Maintain clean work areas.
- Wash hands before leaving the lab.

Instrument Safety

Please refer to the latest Standards for the Safe User of Lasers document from the American National Standards Institute (ANSI Z136.1) for more detailed information to protect against hazards associated with the use of lasers.



WARNING Do not modify this device. Unauthorized modifications may create a safety hazard.










WARNING If this equipment is used in a manner not specified by Fluidigm Corporation, the protection provided by the equipment may be compromised.



WARNING BIOHAZARD. If you are putting biohazardous material on the instrument, use appropriate personal protective equipment and adhere to *Biosafety in Microbiological and Biomedical Laboratories* (BMBL) from the Centers for Disease Control and Prevention and to your lab's safety protocol to limit biohazard risks. If biohazardous materials are used, properly label the equipment as a biohazard. For more information, see the BMBL guidelines online at: [cdc.gov/biosafety/publications/index.htm](https://www.cdc.gov/biosafety/publications/index.htm)

Instrument Symbols

Symbol	English	Description
	Warning	Indicates a hazardous situation that and result in death or serious injury if not avoided.
	Pinch hazard	Indicates where pinch hazards exist. Exercise caution when operating around these areas.
	Laser radiation	The Class 1 laser device is eye-safe under all operating conditions. This product includes lasers of a higher class whose beams are confined within a suitable enclosure so that access to laser radiation is physically prevented.
	Finger Cut Hazard	Broken glass may cause injury or cutting of fingers or hands. Caution when loading and unloading the sample slides.
	Tipping Hazard	Movement or impact with the instrument may cause tipping.
	Trip Hazard	Indicates a tripping hazard. Trip hazards are lying about. Watch your step to avoid falling over objects.
	Compressed gas hazard	Any material contained under pressure, including compressed gas, dissolved gas, or gas liquefied by compression or refrigeration. A compressed gas cylinder can become a projectile when ruptured with the potential to cause significant damage.

Chemical Safety



WARNING Read and comprehend all safety data sheets (SDSs) by chemical manufacturers before you use, store, or handle any chemicals or hazardous materials.



WARNING Wear personal protective equipment (gloves, safety glasses, fully enclosed shoes, lab coats) when handling chemicals.



WARNING Do not inhale fumes from chemicals. Use adequate ventilation, and return caps to bottles immediately after use.



WARNING Check regularly for chemical spills or leaks. Follow SDS recommendations for cleaning up spills or leaks.



For technical support visit fluidigm.com/support