

# Pipetting Maps for the C1 mRNA Seq HT Protocol

## Before You Begin

Before using the maps and for detailed instructions on preparing the C1™ mRNA Seq HT chemistry, see the Using the C1 High-Throughput IFC to Generate Single-Cell cDNA Libraries for mRNA Sequencing Protocol (PN 100-9886).

## Workflow

1 Prime the IFC	2 Load cells	3 Run lysis, reverse transcription, and PCR	4 Harvest cDNA
Prepare pre-mixes, cells, and stain	Pipet and image cells	Prepare mRNA Seq HT chemistry	Transfer cDNA to a 96-well plate for library preparation

## How to Use the Maps

- 1 Select the appropriate map for your workflow step, then print on a color printer as a 1-sided page at actual size (100% scale).
- 2 (Optional) Cut the image at the dotted line (✂️).
- 3 Place the map under the HT IFC to use as a pipetting guide.

## Best Practices

When pipetting into the HT IFC:

- Follow the instructions for each map.
- Ensure that the notch (A1 position) is at the top-left corner of the HT IFC and the barcode faces to the left.
- Always stop at the first stop on the pipette to avoid creating bubbles in the inlets. If a bubble is introduced, ensure that it floats to the top of the well.
- Make sure to evenly distribute reagents over the bottom surfaces of the accumulators and reservoirs.

## Safety

Use standard laboratory safety protocols. Read and understand the safety data sheets (SDSs) before handling chemicals. To obtain SDSs, go to [fluidigm.com/sds](http://fluidigm.com/sds) and search for the SDS using either the product name or the part number.

## For technical support visit [fluidigm.com/support](http://fluidigm.com/support).

**North America** +1 650 266 6100 | Toll-free: 866 358 4354 in the US | [techsupport@fluidigm.com](mailto:techsupport@fluidigm.com)

**Europe** +33 1 60 92 42 40 | [techsupporteurope@fluidigm.com](mailto:techsupporteurope@fluidigm.com)

**China (excluding Hong Kong)** +86 21 3255 8368 | [techsupportchina@fluidigm.com](mailto:techsupportchina@fluidigm.com)

**Japan** +81 3 3662 2150 | [techsupportjapan@fluidigm.com](mailto:techsupportjapan@fluidigm.com)

**All other Asian countries** +1 650 266 6100 | [techsupportasia@fluidigm.com](mailto:techsupportasia@fluidigm.com)






**Central and South America** +1 650 266 6100 | [techsupportlatam@fluidigm.com](mailto:techsupportlatam@fluidigm.com)







## For Research Use Only. Not for use in diagnostic procedures.


Information in this publication is subject to change without notice. **Safety data sheet information:** [fluidigm.com/sds](http://fluidigm.com/sds). **Patent and license information:** [fluidigm.com/legalnotices](http://fluidigm.com/legalnotices). **EU's WEEE directive information:** [fluidigm.com/compliance](http://fluidigm.com/compliance). Fluidigm, the Fluidigm logo, and C1 are trademarks or registered trademarks of Fluidigm Corporation in the United States and/or other countries. All other trademarks are the sole property of their respective owners. © 2016 Fluidigm Corporation. All rights reserved. 01/2016

# Prime the C1 HT IFC

## Reagent Loading

Reagent	Vol. (μL)	Notes
 Valve Fluid	250	The Valve Fluid is a high-surface-tension solution. With a P200 pipette, use 2 aliquots of 125 μL. <b>IMPORTANT</b> Make sure to keep the HT IFC as level as possible so the fluid does not move away from the center of either accumulator (  ). Gently tap the HT IFC on the bench to distribute the fluid evenly over the bottom surface of both accumulators (  ). Use just enough force to <b>avoid creating bubbles</b> .
 Valve Fluid	20	
 Stability Solution	20	The Stability Solution is viscous. The final volume transferred to each well may be <20 μL. <b>IMPORTANT</b> Make sure to use a P200 pipette. Pipet the solution slowly and carefully to <b>avoid bubbles</b> .

Reagent	Vol. (μL)	Notes
 Preloading Reagent	5	Aliquot 30 μL into each tube of an 8-tube strip, then use a multichannel pipette to pipet into each well.
 Preloading Reagent	20	
 Preloading Reagent	150	
 C1 Harvest Reagent	180	
 1X Blocking Reagent	180	<b>IMPORTANT</b> Make sure to dilute the 10X Blocking Reagent to 1X with Cell Rinsing Reagent and keep on ice until use.
 1X Blocking Reagent	20	




# Load Cells into the C1 HT IFC

**IMPORTANT** Remove any remaining reagents in step 1 before loading the reagents in step 2.

## Step 1. Reagent Removal

## Step 2. Reagent Loading

Old Reagent	Previous Map	New Reagent	Vol. (μL)	Notes
■ 1X Blocking Reagent	Prime	□ Cell Rinsing Reagent	180	
● 1X Blocking Reagent	Prime	● • <b>Staining cells:</b> LIVE/DEAD® Staining Solution • <b>Not staining cells:</b> Cell Rinsing Reagent	20	Load one of these reagents.
● 1X Blocking Reagent	Prime	● Cell mix	20	You can load the same or two different cell mixes. <b>IMPORTANT</b> Pipet the cell mix up and down 5–10 times to mix. <b>Do not vortex. Avoid creating bubbles.</b>





# Pipet Lysis, Reverse Transcription, and Preamplification Chemistry into the HT IFC

**IMPORTANT** Remove any remaining reagents in step 1 before loading the reagents in step 2.

## Step 1. Reagent Removal


## Step 2. Reagent Loading

Old Reagent	Previous Map	New Reagent	Vol.( $\mu$ L)	Notes
● Preloading Reagent	Prime	 Lysis Mix A Plus Diluted Barcodes	4.2	Use a multichannel pipette with fresh pipette tips for each column. <b>Avoid creating bubbles.</b> <b>IMPORTANT</b> Do not allow a pipette tip to touch another well.
○ Preloading Reagent	Prime	● RT Mix B	20	
■ Preloading Reagent	Prime	■ Preamplification Mix I (Mix D)	130	
■ C1 Harvest Reagent	Prime	■ Cleanup Mix C	180	If you program to harvest the next day, you can store remaining Cleanup Mix C overnight on ice or at 4 °C until you are ready to harvest.
□ Cell Rinsing Reagent	Load	■ Preloading Reagent	150	



# Harvest cDNA from the C1 HT IFC

## Reagent Removal

Reagent	Vol. (μL)	Notes
 Harvested cDNA amplicons	5	<p>Carefully pull back the tape covering the harvesting inlets on the left and right sides of the HT IFC using the plastic removal tool.</p> <p>Using a multichannel pipette set to 5 μL, transfer the harvested cDNA amplicons from each column of HT IFC inlets and into the corresponding wells of the 96-well harvest plate. Place the plate on ice.</p> <p>The volume harvested from each well may be &lt;5 μL.</p>

